



FSU Biohazardous

Non-Infectious Waste

Non-pathogenic bacteria



Place in **clear, non-labeled** autoclave bag



Waste that is not ready to be autoclaved, ensure that it is in a closed container labeled **"General waste to be autoclaved only"**



Place the clear bag in an autoclaveable container and autoclave on appropriate settings

Sharps

Broken/Unbroken Glass, Pipette Tips, Pasteur Pipettes

Razor Blades, Syringes, Scalpels, Needles



Uncontaminated place in cardboard box labeled **"Broken Glass"** along with building and room number.



Contaminated, place in puncture resistant labeled container



Place in red puncture resistant container



Can be discarded as **Trash**
Should not be dispose as recyclable!



Put in Biomedical Waste Box, discard as **Biomedical Waste** (Not to exceed 30 days)

Infectious Waste

Disease causing organisms, human pathogen containing waste



Place a **red labeled** biohazard autoclave bag



Autoclave on appropriate settings